

## Anti-Zebrafish SUZ12a/b Antibody Cy3 Conjugated

Catalog Number: AZQ6DC03-Cy3

### About SUZ12a/b

Predicted to enable chromatin DNA binding activity. Predicted to act upstream of or within chromatin organization. Predicted to be part of ESC/E(Z) complex and RSC-type complex. Predicted to be active in nucleus. Is expressed in gonad; midbrain hindbrain boundary; pharyngeal arch; and visual system. Human ortholog(s) of this gene implicated in renal Wilms' tumor. Orthologous to human SUZ12 (SUZ12 polycomb repressive complex 2 subunit).

### Overview

Product Name	Anti-Zebrafish SUZ12a/b Antibody Cy3 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6DC03/B0R1D5

### Technical Details

Immunogen	E.coli-derived zebrafish SUZ12a/b recombinant protein (Position: K38-K448).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SUZ12a/b Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.