

Anti-Zebrafish SOX1a/b Antibody Picoband® Cy3 Conjugated

Catalog Number: AZQ6DGL6-Cy3

About SOX1a/b

Contributes to DNA-binding transcription factor activity. Acts upstream of or within determination of left/right symmetry and positive regulation of DNA-binding transcription factor activity. Predicted to be part of transcription regulator complex. Predicted to be active in nucleus. Is expressed in several structures, including forebrain neural keel; lens placode; nervous system; neural plate; and solid lens vesicle. Orthologous to human SOX1 (SRY-box transcription factor 1).

Overview

Product Name	Anti-Zebrafish SOX1a/b Antibody Picoband® Cy3 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6DGL6/Q2Z1R2

Technical Details

Immunogen	E.coli-derived zebrafish SOX1a/b recombinant protein (Position: M185-Y226).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SOX1a/b Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.