

## Anti-Zebrafish SOX1a/b Antibody Picoband® Fluoro550 Conjugated

Catalog Number: AZQ6DGL6-Fluoro550

### About SOX1a/b

Contributes to DNA-binding transcription factor activity. Acts upstream of or within determination of left/right symmetry and positive regulation of DNA-binding transcription factor activity. Predicted to be part of transcription regulator complex. Predicted to be active in nucleus. Is expressed in several structures, including forebrain neural keel; lens placode; nervous system; neural plate; and solid lens vesicle. Orthologous to human SOX1 (SRY-box transcription factor 1).

### Overview

Product Name	Anti-Zebrafish SOX1a/b Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6DGL6/Q2Z1R2

### Technical Details

Immunogen	E.coli-derived zebrafish SOX1a/b recombinant protein (Position: M185-Y226).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SOX1a/b Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.