

Anti-Zebrafish SOX3 Antibody Picoband® HRP Conjugated

Catalog Number: AZQ6EJB7-HRP

About SOX3

Enables DNA-binding transcription activator activity, RNA polymerase II-specific. Contributes to DNA-binding transcription factor activity. Involved in regulation of transcription by RNA polymerase II. Acts upstream of or within several processes, including nervous system development; otic placode development; and otic placode formation. Predicted to be active in nucleus. Is expressed in several structures, including ectoderm; nervous system; neural keel; neural tube; and presumptive structure. Human ortholog(s) of this gene implicated in X-linked panhypopituitarism. Orthologous to human SOX3 (SRY-box transcription factor 3).

Overview

Product Name	Anti-Zebrafish SOX3 Antibody Picoband® HRP Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6EJB7

Technical Details

Immunogen	E.coli-derived Zebrafish SOX3 recombinant protein (Position: Q141-Q285)
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SOX3 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.