

Anti-Zebrafish FABP6 Antibody Picoband® FITC Conjugated

Catalog Number: AZQ6IMW5-FITC

About FABP6

Enables bile acid binding activity. Predicted to be involved in fatty acid transport. Predicted to be located in cytoplasm and membrane. Predicted to be active in cytosol and nucleus. Is expressed in digestive system; heart; and pleuroperitoneal region. Orthologous to human FABP6 (fatty acid binding protein 6).

Overview

Product Name	Anti-Zebrafish FABP6 Antibody Picoband® FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6IMW5

Technical Details

Immunogen	E.coli-derived Zebrafish FABP6 recombinant protein (Position: P25-V131).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish FABP6 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.