

Anti-Zebrafish SNRPD3/3l Antibody Picoband® FITC Conjugated

Catalog Number: AZQ6IQ56-FITC

About SNRPD3/3l

Predicted to enable RNA binding activity. Predicted to be involved in spliceosomal snRNP assembly. Predicted to act upstream of or within RNA splicing and mRNA processing. Predicted to be located in cytosol and nucleus. Predicted to be part of SMN-Sm protein complex; spliceosomal complex; and spliceosomal snRNP complex. Human ortholog(s) of this gene implicated in systemic lupus erythematosus. Orthologous to human SNRPD3 (small nuclear ribonucleoprotein D3 polypeptide).

Overview

Product Name	Anti-Zebrafish SNRPD3/3l Antibody Picoband® FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6IQ56/Q7ZVB5

Technical Details

Immunogen	E.coli-derived zebrafish SNRPD3/3l recombinant protein (Position: M1-R112).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SNRPD3/3I Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.