

## Anti-Zebrafish TSG101a/b Antibody Picoband® Fluoro488 Conjugated

Catalog Number: AZQ6IQ70-Fluoro488

### About TSG101a/b

Predicted to enable ubiquitin binding activity. Predicted to be involved in endosome to lysosome transport. Predicted to act upstream of or within protein modification process. Predicted to be located in endosome. Predicted to be part of ESCRT I complex. Human ortholog(s) of this gene implicated in breast cancer; breast carcinoma; cervix carcinoma; endometrial adenocarcinoma; and prostate cancer. Orthologous to human TSG101 (tumor susceptibility 101).

### Overview

Product Name	Anti-Zebrafish TSG101a/b Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6IQ70/E7F4Y8

### Technical Details

Immunogen	E.coli-derived zebrafish TSG101a/b recombinant protein (Position: G222-S275).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



**Anti-Zebrafish TSG101a/b Antibody - Fluoro488**

For Research Use Only. Not for use in diagnostic procedures.