

Anti-Zebrafish GATA6 Antibody Biotin Conjugated

Catalog Number: AZQ6NW63-Biotin

About GATA6

Enables DNA binding activity. Acts upstream of or within several processes, including cardioblast migration to the midline involved in heart field formation; liver development; and positive regulation of transcription by RNA polymerase II. Predicted to be active in nucleus. Is expressed in several structures, including blastoderm; cardiovascular system; digestive system; mesoderm; and yolk syncytial layer. Human ortholog(s) of this gene implicated in adenocarcinoma (multiple); congenital heart disease (multiple); and pancreatic hypoplasia-diabetes-congenital heart disease syndrome. Orthologous to human GATA6 (GATA binding protein 6).

Overview

Product Name	Anti-Zebrafish GATA6 Antibody Biotin Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6NW63

Technical Details

Immunogen	E.coli-derived Zebrafish GATA6 recombinant protein (Position: E26-E296)
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish GATA6 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.