

Anti-Zebrafish SF3B4 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: AZQ6NWB3-Fluoro550

About SF3B4

Predicted to enable RNA binding activity. Acts upstream of or within regulation of retina development in camera-type eye and visual behavior. Predicted to be located in nucleus. Predicted to be part of U2 snRNP and precatalytic spliceosome. Is expressed in head. Human ortholog(s) of this gene implicated in Nager acrofacial dysostosis; cervical cancer; and hepatocellular carcinoma. Orthologous to human SF3B4 (splicing factor 3b subunit 4).

Overview

Product Name	Anti-Zebrafish SF3B4 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NWB3

Technical Details

Immunogen	E.coli-derived zebrafish SF3B4 recombinant protein (Position: M1-A214).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SF3B4 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.