

Anti-Zebrafish Cytokeratin 8/KRT8 Antibody Picoband® APC Conjugated

Catalog Number: AZQ6NWF6-APC

About KRT8

Predicted to be a structural constituent of skin epidermis. Acts upstream of or within camera-type eye development and neural crest cell migration. Predicted to be located in cytoplasm. Predicted to be active in keratin filament. Is expressed in several structures, including EVL; cardiovascular system; integument; pectoral fin; and sensory system. Human ortholog(s) of this gene implicated in several diseases, including Meesmann corneal dystrophy 2; autosomal dominant woolly hair; hypotrichosis (multiple); ichthyosis (multiple); and palmoplantar keratosis (multiple). Orthologous to several human genes including KRT8 (keratin 8).

Overview

Product Name	Anti-Zebrafish Cytokeratin 8/KRT8 Antibody Picoband® APC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NWF6

Technical Details

Immunogen	E.coli-derived Zebrafish Cytokeratin 8/KRT8 recombinant protein (Position: S17-D520).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Zebrafish Cytokeratin 8/KRT8 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.