

## Anti-Zebrafish RCC2 Antibody Picoband® FITC Conjugated

Catalog Number: AZQ6NYE2-FITC

### About RCC2

Predicted to enable small GTPase binding activity. Involved in neural crest cell migration. Predicted to be located in several cellular components, including chromosome, centromeric region; midbody; and spindle. Predicted to be active in membrane. Orthologous to human RCC2 (regulator of chromosome condensation 2).

### Overview

Product Name	Anti-Zebrafish RCC2 Antibody Picoband® FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NYE2

### Technical Details

Immunogen	E.coli-derived zebrafish RCC2 recombinant protein (Position: K266-L495).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish RCC2 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.