

Anti-Zebrafish RCC2 Antibody Picoband® HRP Conjugated

Catalog Number: AZQ6NYE2-HRP

About RCC2

Predicted to enable small GTPase binding activity. Involved in neural crest cell migration. Predicted to be located in several cellular components, including chromosome, centromeric region; midbody; and spindle. Predicted to be active in membrane. Orthologous to human RCC2 (regulator of chromosome condensation 2).

Overview

Product Name	Anti-Zebrafish RCC2 Antibody Picoband® HRP Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6NYE2

Technical Details

Immunogen	E.coli-derived zebrafish RCC2 recombinant protein (Position: K266-L495).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.