

Anti-Zebrafish TDRD3 Antibody Picoband® FITC Conjugated

Catalog Number: AZQ6NYG6-FITC

About TDRD3

Predicted to enable chromatin binding activity; methylated histone binding activity; and transcription coactivator activity. Predicted to act upstream of or within chromatin organization and mRNA processing. Predicted to be located in cytoplasm. Predicted to be part of exon-exon junction complex. Predicted to be active in nucleus. Orthologous to human TDRD3 (tudor domain containing 3).

Overview

Product Name	Anti-Zebrafish TDRD3 Antibody Picoband® FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NYG6

Technical Details

Immunogen	E.coli-derived zebrafish TDRD3 recombinant protein (Position: K216-R242).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish TDRD3 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.