

## Anti-Zebrafish TDRD3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: AZQ6NYG6-Fluoro488

### About TDRD3

Predicted to enable chromatin binding activity; methylated histone binding activity; and transcription coactivator activity. Predicted to act upstream of or within chromatin organization and mRNA processing. Predicted to be located in cytoplasm. Predicted to be part of exon-exon junction complex. Predicted to be active in nucleus. Orthologous to human TDRD3 (tudor domain containing 3).

### Overview

Product Name	Anti-Zebrafish TDRD3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NYG6

### Technical Details

Immunogen	E.coli-derived zebrafish TDRD3 recombinant protein (Position: K216-R242).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



**Anti-Zebrafish TDRD3 Antibody - Fluoro488**

For Research Use Only. Not for use in diagnostic procedures.