

Anti-Zebrafish CUL1A/B Antibody Picoband® Fluoro647 Conjugated

Catalog Number: AZQ6NYS3-Fluoro647

Overview

Product Name	Anti-Zebrafish CUL1A/B Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NYS3,Q802D4

Technical Details

Immunogen	E.coli-derived zebrafish CUL1A/B recombinant protein (Position: K436-A777).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish CUL1A/B Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.