

Anti-Zebrafish LEO1 Antibody Picoband® FITC Conjugated

Catalog Number: AZQ6NYV9-FITC

About LEO1

RNA polymerase-associated protein LEO1 is an enzyme that in humans is encoded by the LEO1 gene. LEO1, parafibromin (CDC73; MIM 607393), CTR9 (MIM 609366), and PAF1 (MIM 610506) form the PAF protein complex that associates with the RNA polymerase II subunit POLR2A (MIM 180660) and with a histone methyltransferase complex

Overview

Product Name	Anti-Zebrafish LEO1 Antibody Picoband® FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6NYV9

Technical Details

Immunogen	E.coli-derived zebrafish LEO1 recombinant protein (Position: E333-K571).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish LEO1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.