

Anti-Zebrafish OSTF1 Antibody Picoband® Biotin Conjugated

Catalog Number: AZQ6TGW5-Biotin

About OSTF1

Osteoclast-stimulating factor 1 is a protein that in humans is encoded by the OSTF1 gene. Osteoclast-stimulating factor-1 is an intracellular protein produced by osteoclasts that indirectly induces osteoclast formation and bone resorption.

Overview

Product Name	Anti-Zebrafish OSTF1 Antibody Picoband® Biotin Conjugated
Reactive Species	Zebrafish
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6TGW5

Technical Details

Immunogen	E.coli-derived zebrafish OSTF1 recombinant protein (Position: R57-K190).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish OSTF1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.