

## Anti-Zebrafish SF3A1 Antibody FITC Conjugated

Catalog Number: AZQ90X41-FITC

### About SF3A1

Predicted to enable RNA binding activity. Predicted to act upstream of or within mRNA cis splicing, via spliceosome. Predicted to be located in nucleus. Predicted to be part of U2 snRNP; U2-type prespliceosome; and catalytic step 2 spliceosome. Orthologous to human SF3A1 (splicing factor 3a subunit 1).

### Overview

Product Name	Anti-Zebrafish SF3A1 Antibody FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (IF, IHC). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q90X41

### Technical Details

Immunogen	E.coli-derived zebrafish SF3A1 recombinant protein (Position: K16-P646).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SF3A1 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.