

Anti-Zebrafish SF3A1 Antibody Fluoro550 Conjugated

Catalog Number: AZQ90X41-Fluoro550

About SF3A1

Predicted to enable RNA binding activity. Predicted to act upstream of or within mRNA cis splicing, via spliceosome. Predicted to be located in nucleus. Predicted to be part of U2 snRNP; U2-type prespliceosome; and catalytic step 2 spliceosome. Orthologous to human SF3A1 (splicing factor 3a subunit 1).

Overview

Product Name	Anti-Zebrafish SF3A1 Antibody Fluoro550 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q90X41

Technical Details

Immunogen	E.coli-derived zebrafish SF3A1 recombinant protein (Position: K16-P646).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish SF3A1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.