

Anti-Zebrafish Cadherin-17/CDH17 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: AZQ90X63-Fluoro594

About CDH17

Predicted to enable beta-catenin binding activity and cadherin binding activity. Acts upstream of or within pronephros development. Predicted to be located in plasma membrane. Predicted to be part of catenin complex. Predicted to be active in adherens junction. Is expressed in digestive system; intermediate cell mass of mesoderm; mesoderm; and renal system. Orthologous to human CDH17 (cadherin 17).

Overview

Product Name	Anti-Zebrafish Cadherin-17/CDH17 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q90X63

Technical Details

Immunogen	E.coli-derived Zebrafish Cadherin-17/CDH17 recombinant protein (Position: E195-L868).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish Cadherin-17/CDH17 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.