

Anti-Zebrafish JAG1B Antibody Picoband® HRP Conjugated

Catalog Number: AZQ90Y54-HRP

About JAG1B

Predicted to enable Notch binding activity. Acts upstream of or within several processes, including exocrine pancreas development; face morphogenesis; and inner ear development. Predicted to be located in membrane. Is expressed in several structures, including digestive system; endocrine system; nervous system; otic vesicle; and pharyngeal arch. Human ortholog(s) of this gene implicated in Alagille syndrome; Charcot-Marie-Tooth disease; congenital heart disease; multiple sclerosis; and tetralogy of Fallot. Orthologous to human JAG1 (jagged canonical Notch ligand 1).

Overview

Product Name	Anti-Zebrafish JAG1B Antibody Picoband® HRP Conjugated
Reactive Species	Zebrafish
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q90Y54

Technical Details

Immunogen	E.coli-derived Zebrafish JAG1B recombinant protein (Position: Q853-V1213).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish JAG1B Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.