

Anti-Zebrafish FZD2 Antibody Picoband® FITC Conjugated

Catalog Number: AZQ90YL7-FITC

About FZD2

Predicted to enable Wnt receptor activity and Wnt-protein binding activity. Acts upstream of or within several processes, including Wnt signaling pathway, calcium modulating pathway; heart development; and sensory epithelium regeneration. Located in plasma membrane. Is expressed in several structures, including cardiovascular system; fin; head; head mesenchyme; and mesoderm. Human ortholog(s) of this gene implicated in omodyspasia 2. Orthologous to human FZD2 (frizzled class receptor 2).

Overview

Product Name	Anti-Zebrafish FZD2 Antibody Picoband® FITC Conjugated
Reactive Species	Zebrafish
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q90YL7

Technical Details

Immunogen	E.coli-derived Zebrafish FZD2 recombinant protein (Position: E134-V550).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish FZD2 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.