

Anti-Zebrafish FZD2 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: AZQ90YL7-Fluoro647

About FZD2

Predicted to enable Wnt receptor activity and Wnt-protein binding activity. Acts upstream of or within several processes, including Wnt signaling pathway, calcium modulating pathway; heart development; and sensory epithelium regeneration. Located in plasma membrane. Is expressed in several structures, including cardiovascular system; fin; head; head mesenchyme; and mesoderm. Human ortholog(s) of this gene implicated in omodysplasia 2. Orthologous to human FZD2 (frizzled class receptor 2).

Overview

Product Name	Anti-Zebrafish FZD2 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q90YL7

Technical Details

Immunogen	E.coli-derived Zebrafish FZD2 recombinant protein (Position: E134-V550).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish FZD2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.