

Anti-Zebrafish Bcl-XL/BCL2L1 Antibody Fluoro488 Conjugated

Catalog Number: AZQ90Z98-Fluoro488

About BCL2L1

Enables cysteine-type endopeptidase inhibitor activity involved in apoptotic process. Acts upstream of or within negative regulation of intrinsic apoptotic signaling pathway; thyroid gland development; and tissue homeostasis. Predicted to be located in membrane. Predicted to be active in mitochondrial outer membrane. Is expressed in several structures, including blood island; endocrine system; extension; immune system; and nervous system. Human ortholog(s) of this gene implicated in breast cancer; carcinoma (multiple); leukemia (multiple); and pancreatic cancer. Orthologous to human BCL2L1 (BCL2 like 1).

Overview

Product Name	Anti-Zebrafish Bcl-XL/BCL2L1 Antibody Fluoro488 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q90Z98

Technical Details

Immunogen	E.coli-derived Zebrafish Bcl-XL/BCL2L1 recombinant protein (Position: Q19-K215).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-Zebrafish Bcl-XL/BCL2L1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.