

Anti-Zebrafish NKX2.1 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: AZQ98TX4-Fluoro647

About NKX2.1

Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Acts upstream of or within hypothalamus development. Predicted to be active in nucleus. Is expressed in brain; head; presumptive diencephalon; swim bladder; and thyroid follicle. Human ortholog(s) of this gene implicated in choreatic disease and papillary thyroid carcinoma. Orthologous to human NKX2-1 (NK2 homeobox 1).

Overview

Product Name	Anti-Zebrafish NKX2.1 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q98TX4

Technical Details

Immunogen	E.coli-derived Zebrafish NKX2.1 recombinant protein (Position: S23-S326).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish NKX2.1 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.