

## Anti-Zebrafish SST2 Antibody Cy3 Conjugated

Catalog Number: AZQ9DDE4-Cy3

### About SST2

Predicted to enable hormone activity. Predicted to be involved in regulation of cell migration. Predicted to be located in extracellular region. Predicted to be active in extracellular space. Is expressed in brainstem; endocrine system; floor plate; and pleuroperitoneal region.

### Overview

Product Name	Anti-Zebrafish SST2 Antibody Cy3 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9DDE4

### Technical Details

Immunogen	E.coli-derived Zebrafish SST2 recombinant protein (Position: R25-C107).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.