

Anti-Zebrafish VENT Antibody Picoband® Fluoro550 Conjugated

Catalog Number: AZQ9DGM9-Fluoro550

About VENT

Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Acts upstream of or within dorsal/ventral pattern formation; regulation of endodermal cell fate specification; and regulation of mesodermal cell fate specification. Predicted to be active in nucleus. Is expressed in embryonic structure. Orthologous to human VENTX (VENT homeobox).

Overview

Product Name	Anti-Zebrafish VENT Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Zebrafish
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9DGM9

Technical Details

Immunogen	E.coli-derived Zebrafish VENT recombinant protein (Position: M1-F170).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Zebrafish VENT Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.