

Anti-Saccharomyces cerevisiae TY1RT239-539 Antibody HRP Conjugated

Catalog Number: DZ33991-HRP

Description

Size	200 ul/vial
Form	Liquid
Isotype	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Saccharomyces cerevisiae TY1RT239-539 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.