

Anti-Saccharomyces cerevisiae Ty1RT217-234 Antibody Biotin Conjugated

Catalog Number: DZ33992-Biotin

Description

Size	200 ul/vial
Form	Liquid
Isotype	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Saccharomyces cerevisiae Ty1RT217-234 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.