

Anti-Red wild einkorn MLH1 TRIUR3_24234 Antibody Biotin Conjugated

Catalog Number: DZ41013-Biotin

Description

Size	200 ul/vial
Form	Liquid
Isotype	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Red wild einkorn MLH1 TRIUR3_24234 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.