

## Anti-Nurse Shark b2m Antibody Fluoro488 Conjugated

Catalog Number: DZ41070-Fluoro488

### Description

|                     |  |
|---------------------|--|
| Size                | 200 ul/vial  |
| Form                | Liquid   |
| Isotype             | Rabbit IgG   |
| Purification        | Immunogen affinity purified.   |
| Contents            | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Conjugate           | Fluoro488<br>Excitation Wavelength: 488 nm<br>Emission Wavelength: 515-545 nm                                |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.   |

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nurse Shark b2m Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.