

Anti-Nurse Shark class II beta Antibody Biotin Conjugated

Catalog Number: DZ41071-Biotin

Description

| | |
|---------------------|---|
| Size | 200 ul/vial |
| Form | Liquid |
| Isotype | Rabbit IgG |
| Purification | Immunogen affinity purified. |
| Contents | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Conjugate | Biotin |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nurse Shark class II beta Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.