

## Anti-African malaria mosquito CTLMA2 Antibody Biotin Conjugated

Catalog Number: DZ41078-Biotin

### Description

|                     |   |
|---------------------|---|
| Size                | 200 ul/vial   |
| Form                | Liquid  |
| Isotype             | Rabbit IgG  |
| Purification        | Immunogen affinity purified.  |
| Contents            | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .  |
| Conjugate           | Biotin  |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users.<br>Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.<br>ELISA, Optimal dilutions should be determined by end users. |

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-African malaria mosquito CTLMA2 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.