

Anti-African malaria mosquito LRRD7 Antibody Fluoro488 Conjugated

Catalog Number: DZ41082-Fluoro488

Description

Size	200 ul/vial
Form	Liquid
Isotype	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-African malaria mosquito LRRD7 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.