

## Anti-African malaria mosquito Rel2 Antibody Fluoro647 Conjugated

Catalog Number: DZ41083-Fluoro647

### Description

Size	200 ul/vial
Form	Liquid
Isotype	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-African malaria mosquito Rel2 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.