

Anti-Yellowfever mosquito OBP22 Antibody Fluoro647 Conjugated

Catalog Number: DZ41084-Fluoro647

Description

| | |
|---------------------|--|
| Size | 200 ul/vial |
| Form | Liquid |
| Isotype | Rabbit IgG |
| Purification | Immunogen affinity purified. |
| Contents | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Yellowfever mosquito OBP22 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.