

## Anti-Austrofundulus limnaeus g3bp1 Antibody HRP Conjugated

Catalog Number: DZ41405-HRP

### Description

Size	200 ul/vial
Form	Liquid
Isotype	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Austrofundulus limnaeus g3bp1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.