

Certificate of Analysis

Product Name: Polyclonal Anti-Drosophila melanogaster (Fruit fly) RE08173p/Dmel/Tg Antibody

Catalog Number: DZ41557

Description

Lot No.	24DZ670F13
Size	200µl/vial
Form	Liquid
Concentration	500µg/ml
Ig type	Rabbit IgG
Purification	Immunogen affinity purified.
Contents	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN ₃ .

Data of Indirect ELISA

Coating concentration of antigen	0.1µg/ml	0.1µg/ml	0.1µg/ml	0.1µg/ml
Working concentration of anti-BRE08173p/Dmel/Tg antibody	1µg/ml	0.1µg/ml	0.01µg/ml	0.001µg/ml
Working concentration of HRP conjugated goat anti-rabbit IgG (Catalog # BA1054)	0.2µg/ml	0.2µg/ml	0.2µg/ml	0.2µg/ml
O.D. values(450nm)	2.792	2.312	0.495	0.137
Antibody Status	Pass			


Data Evaluation Criteria

Coating concentration of antigen	0.1µg/ml	0.1µg/ml
Working concentration of the primary antibody	1µg/ml	1µg/ml
Working concentration of HRP conjugated goat anti-rabbit IgG	0.2µg/ml	0.2µg/ml
O.D. values (450nm)	≥ 1.0	<1.0
Antibody's affinity for the antigen	Strong	Weak
Antibody Status	Pass	Fail

Storage

At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.

Approved by:



Ms. Yangjing Yue

Chief technician

Quality Control System of ELISA

Date: June 13, 2024