

Anti-Tribolium castaneum (Red flour beetle) AUGUSTUS-3.0.2_09181 Antibody Fluoro647 Conjugated

Catalog Number: DZ41686-Fluoro647

Description

Size	200 ul/vial
Form	Liquid
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tribolium castaneum (Red flour beetle) AUGUSTUS-3.0.2_09181 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.