

Anti-Parhyale hawaiiensis BMAL1 Antibody

Catalog Number: DZ41687

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Parhyale hawaiiensis BMAL1 Antibody
Reactive Species	Parhyale hawaiiensis
Description	Boster Bio Anti-Parhyale hawaiiensis BMAL1 Antibody catalog # DZ41687. Tested in WB, IHC applications. This antibody reacts with Parhyale hawaiiensis.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Host	Rabbit

Technical Details

Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml Immunohistochemistry(Paraffin-embedded Section), 2-5 ug/ml

[Anti-Parhyale hawaiiensis BMAL1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.