

## Anti-Parhyale hawaiiensis PER Antibody Fluoro594 Conjugated

Catalog Number: DZ41689-Fluoro594

### Description

|                     |  |
|---------------------|--|
| Size                | 200 ul/vial  |
| Form                | Liquid   |
| Purification        | Immunogen affinity purified.   |
| Contents            | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N. |
| Conjugate           | Fluoro594<br>Excitation Wavelength: 593 nm<br>Emission Wavelength: 618 nm                                    |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.   |

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Parhyale hawaiiensis PER Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.