

## Anti-Tomato yellow leaf curl virus (TYLCV) Protein C4 Antibody

Catalog Number: DZ41692

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview

|                      |                                                                                                                                                              |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Tomato yellow leaf curl virus (TYLCV) Protein C4 Antibody                                                                                               |
| Reactive Species     | TYLCV                                                                                                                                                        |
| Description          | Boster Bio Anti-Tomato yellow leaf curl virus (TYLCV) Protein C4 Antibody catalog # DZ41692. Tested in WB, IF applications. This antibody reacts with TYLCV. |
| Application          | IF, WB                                                                                                                                                       |
| Clonality            | Polyclonal                                                                                                                                                   |
| Formulation          | Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN <sub>3</sub> .                                                                                          |
| Storage Instructions | At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.                                                                            |
| Host                 | Rabbit                                                                                                                                                       |
| Uniprot ID           | P27271                                                                                                                                                       |

### Technical Details

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Form                | Liquid                                                      |
| Concentration       | 500 ug/ml                                                   |
| Purification        | Immunogen affinity purified.                                |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml<br>Immunofluorescence, 5 ug/ml |

Anti-Tomato yellow leaf curl virus (TYLCV) Protein C4 Antibody

For Research Use Only. Not for use in diagnostic procedures.