

## Anti-Chlamydomonas reinhardtii (Chlamydomonas smithii) plus agglutinin/SAG1 Antibody Fluoro594 Conjugated

Catalog Number: DZ41707-Fluoro594

### Description

Size	200 ul/vial
Form	Liquid
Purification	Immunogen affinity purified.
Contents	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Chlamydomonas reinhardtii (Chlamydomonas smithii) plus agglutinin/SAG1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.