

Certificate of Analysis

Product Name: Polyclonal Anti-Chlamydomonas reinhardtii (Chlamydomonas smithii)
plus agglutinin/SAG1 Antibody

Catalog Number: DZ41707

Description

| | |
|----------------------|---|
| Lot No. | 25DZ835D29 |
| Size | 200µl/vial |
| Form | Liquid |
| Concentration | 500µg/ml |
| Ig type | Rabbit IgG |
| Purification | Immunogen affinity purified. |
| Contents | Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN ₃ . |

Data of Indirect ELISA

| | | | | |
|---|----------|----------|-----------|------------|
| Coating concentration of antigen | 0.1µg/ml | 0.1µg/ml | 0.1µg/ml | 0.1µg/ml |
| Working concentration of anti-plus agglutinin/SAG1 antibody | 1µg/ml | 0.1µg/ml | 0.01µg/ml | 0.001µg/ml |
| Working concentration of HRP conjugated goat anti-rabbit IgG (Catalog # BA1054) | 0.2µg/ml | 0.2µg/ml | 0.2µg/ml | 0.2µg/ml |
| O.D. values(450nm) | 4.8060 | 4.1687 | 2.4178 | 0.7460 |
| Antibody Status | Pass | | | |


Data Evaluation Criteria

| | | |
|--|----------|----------|
| Coating concentration of antigen | 0.1µg/ml | 0.1µg/ml |
| Working concentration of the primary antibody | 1µg/ml | 1µg/ml |
| Working concentration of HRP conjugated goat anti-rabbit IgG | 0.2µg/ml | 0.2µg/ml |
| O.D. values (450nm) | ≥ 1.0 | <1.0 |
| Antibody's affinity for the antigen | Strong | Weak |
| Antibody Status | Pass | Fail |

Storage

At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.

Approved by:



Ms. Yangjing Yue

Chief technician

Quality Control System of ELISA

Date: Apr 29, 2025