

Anti-Citrus fleck mosaic virus capsid protein Antibody HRP Conjugated

Catalog Number: DZ41788-HRP

Description

| | |
|---------------------|---|
| Size | 200 ul/vial |
| Form | Liquid |
| Purification | Immunogen affinity purified. |
| Contents | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Conjugate | HRP |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users. |

Storage

At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Citrus fleck mosaic virus capsid protein Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.