

Anti-Citrus fleck mosaic virus capsid protein Antibody

Catalog Number: DZ41788-carrier-free

Overview

Product Name	Anti-Citrus fleck mosaic virus capsid protein Antibody
Reactive Species	Citrus fleck mosaic virus
Description	Boster Bio Anti-Citrus fleck mosaic virus capsid protein Antibody catalog # DZ41788. This antibody reacts with Citrus fleck mosaic virus.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN3.
Storage Instructions	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Host	Rabbit

Technical Details

Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml ELISA, 0.1-0.5 ug/ml

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Citrus fleck mosaic virus capsid protein Antibody

For Research Use Only. Not for use in diagnostic procedures.