

Anti- Aedes aegypti (Yellowfever mosquito) AeLoqs Antibody

Catalog Number: DZ41822

Overview

| | |
|----------------------|--|
| Product Name | Anti- Aedes aegypti (Yellowfever mosquito) AeLoqs Antibody |
| Reactive Species | Yellowfever mosquito |
| Description | Boster Bio Anti- Aedes aegypti (Yellowfever mosquito) AeLoqs Antibody catalog # DZ41822. This antibody reacts with Yellowfever mosquito. |
| Application | IF, IHC, WB |
| Clonality | Polyclonal |
| Formulation | Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN3. |
| Storage Instructions | At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | A0A0D3QFR7 |

Technical Details

| | |
|---------------------|--|
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml Immunohistochemistry (Paraffin-embedded Section), 2-5 ug/ml Immunofluorescence, 5 ug/ml |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti- Aedes aegypti (Yellowfever mosquito) AeLoqs Antibody

For Research Use Only. Not for use in diagnostic procedures.