

## Anti-Aedes aegypti (Yellowfever mosquito) AeDc2 Antibody

Catalog Number: DZ41824

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-Aedes aegypti (Yellowfever mosquito) AeDc2 Antibody   |
| Reactive Species     | Yellowfever mosquito   |
| Description          | Boster Bio Anti-Aedes aegypti (Yellowfever mosquito) AeDc2 Antibody catalog # DZ41824. This antibody reacts with Yellowfever mosquito. |
| Application          | IF, IHC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN3.   |
| Storage Instructions | At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.  |
| Host                 | Rabbit   |
| Uniprot ID           | A0A1S4FEL6   |

### Technical Details

|                     |  |
|---------------------|--|
| Form                | Liquid   |
| Concentration       | 500 ug/ml  |
| Purification        | Immunogen affinity purified.   |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml<br>Immunohistochemistry (Paraffin-embedded Section), 2-5 ug/ml<br>Immunofluorescence, 5 ug/ml |

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Aedes aegypti (Yellowfever mosquito) AeDc2 Antibody

For Research Use Only. Not for use in diagnostic procedures.