

## Anti-Anopheles gambiae (African malaria mosquito) AeA0A6I8T3J6 Antibody

Catalog Number: DZ41830

### Overview

Product Name	Anti-Anopheles gambiae (African malaria mosquito) AeA0A6I8T3J6 Antibody
Reactive Species	African malaria mosquito
Description	Boster Bio Anti-Anopheles gambiae (African malaria mosquito) AeA0A6I8T3J6 Antibody catalog # DZ41830. This antibody reacts with African malaria mosquito.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	A0A6I8T3J6

### Technical Details

Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml Immunohistochemistry (Paraffin-embedded Section), 2-5 ug/ml Immunofluorescence, 5 ug/ml

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Anopheles gambiae (African malaria mosquito) AeA0A6I8T3J6 Antibody

For Research Use Only. Not for use in diagnostic procedures.