

Anti-Anopheles gambiae (African malaria mosquito) AgAgMDL8 Antibody

Catalog Number: DZ41842

Overview

| | |
|----------------------|---|
| Product Name | Anti-Anopheles gambiae (African malaria mosquito) AgAgMDL8 Antibody |
| Reactive Species | African malaria mosquito |
| Description | Boster Bio Anti-Anopheles gambiae (African malaria mosquito) AgAgMDL8 Antibody catalog # DZ41842. This antibody reacts with African malaria mosquito. |
| Application | ELISA, IF, WB |
| Clonality | Polyclonal |
| Formulation | Each vial contains 20mM PBS, 50% glycerol, 0.02% NaN3. |
| Storage Instructions | At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing. |
| Host | Rabbit |

Technical Details

| | |
|---------------------|---|
| Form | Liquid |
| Concentration | 500 ug/ml |
| Purification | Immunogen affinity purified. |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml Immunofluorescence, 5 ug/ml ELISA, 0.1-0.5 ug/ml |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Anopheles gambiae (African malaria mosquito) AgAgMDL8 Antibody

For Research Use Only. Not for use in diagnostic procedures.